

Product datasheet

Baculo Expression Vector pAc-lambda-Fc

Short overview

Cat. No.	PR3003
Quantity	5 µg
Concentration	0.5 µg/µl in TE buffer

Product description

Host	E. coli
Purity	OD 260/280 ratio: 1.8-2.0
Quality check	Efficiency in standard transformation procedure with chemically competent E. coli cells (XL1-Blue) > 1E+08 cfu/µg DNA
Application	Cloning of heavy and light chain variable domain genes of scFv antibody fragments selected from hybridomas, single B-lymphocytes or phage display libraries; Expression of immunoglobulin heavy and light chains in insect cells and secretion of assembled IgG(kappa) or IgG(lamba) antibodies into the supernatant
Stability	Minimum 1 year when stored at -20°C
Vector	pAcUW51 derivative, AmpR, 6701 bp
Cloning site	XhoI and SpeI for heavy chain Fab fd genes, SacI and EcoRV for light chain Fab (lambda) genes
Purification	Plasmid Purification Kit
Storage	-20°C
Intended use	Research use only

Background

Introduction Recombinant antibody technologies have been widely used to produce various single-chain Fv or Fab antibody fragments of different specificity. The randomized combination of cloned variable heavy and light chain immunoglobulin gene fragments further allowed the construction of human antibody libraries, which enable today the isolation of specific scFv or Fab antibody fragments against particular antigens e.g. by phage display. For many applications, however, it is required to reassemble the variable regions of the selected antibody with immunoglobulin constant regions to generate complete antibody molecules. The baculovirus expression system has been established as a reliable system for the production of immunoglobulins (Nesbit et al., 1992; Liang et al., 1997, 2001). Even the development of cassette vectors for the production of human antibodies have been reported (Poul et al., 1995a,b). However, these systems are based on a combination of two different vectors which separately served for the expression of light and heavy chains. The required careful and time consuming adjustment of the two respective recombinant baculovirus titers becomes obsolete with the use of PROGEN's single baculovirus cassette vectors with authentic IgG heavy and light chain signal sequences for the rapid production of complete chimeric, humanized and human IgG antibodies in recombinant baculovirus infected insect cells (Liang et al., 2001). The vectors were specifically designed for cloning of heavy and light chain Fv or Fab gene fragments isolated from hybridomas, individual B cell clones as well as antibody libraries. Vector Design To achieve the expression of complete immunoglobulins, PROGEN offers a set of baculovirus cassette vectors for the convenient insertion of heavy and light chain Fv or Fab domain coding regions. The vectors are based on the backbone of pAcUW51 (BD Biosciences), which contains the f1 origin of phage DNA replication, an ampicillin resistance gene for selection in E. coli, the two baculovirus expression promoters of polyhedrin and p10 and SV40 transcription terminators. The gene elements for IgG expression were cloned successively into pAcUW51. The heavy chain gene cassette preceded by the authentic IgG signal sequence from IgG1 (subgroup VHIII) is located under control of the polyhedrin promoter. In pAc-kappa-CH3 and

