

Product datasheet

anti-Yersinia enterocolitica O:9 mouse monoclonal, 8E9, FITC Conjugate

Short overview

Cat. No.	8E9-FITC
Quantity	1 ml

Product description

Mouse
Monoclonal
IgG3
8E9
Yersinia enterocolitica serogroup O:9 strain
Contains 10% FBS and 0.09% sodium azide in PBS
FITC
Affinity chromatography
2-8°C
Research use only
ICC/IF
Y. enterocolitica O:9

Applications


Immunocytochemistry (ICC)

1:20


Background

8E9 specifically reacts with Yersinia enterocolitica serogroup O:9 strains in both agglutination and immunofluorescence assays. The antibody does not react with Y. enterocolitica O:3 and O:8 serogroups.

Product images



Immunofluorescence analysis of Yersinia enterocolitica infected EMT6 mouse breast cancer cells with anti-Yersinia enterocolitica O:9. Cell were infected with defined amounts of Yersinia enterocolitica (multiplicity of infection MOI 50) for 45 minutes. Fixation was performed using either 4% paraformaldehyde (PFA), 1% PFA, 100% methanol or 100% acetone for 10 min at RT. Cells were blocked and permeabilized with 1% BSA and 0.3% Triton-X100 in PBS for 1 h at RT. The primary antibody anti-Yersinia enterocolitica O:9 mouse monoclonal, 8E9, FITC Conjugate (Cat. No. 8E9-FITC) was diluted in blocking buffer (1:20) and incubated over-night at 4°C. DNA was stained with Hoechst in blue.



Immunofluorescence analysis of Yersinia enterocolitica infected EMT6 mouse breast cancer cells with anti-Yersinia enterocolitica O:9. Cell were infected with defined amounts of Yersinia enterocolitica (multiplicity of infection MOI 50, MOI 5 or without bacteria) for 45 minutes. Fixation was performed using 1% paraformaldehyde for 10 min at RT. Cells were blocked and permeabilized with 1% BSA and 0.3% Triton-X100 in PBS for 1 h at RT. The primary antibody anti-Yersinia enterocolitica O:9 mouse monoclonal, 8E9, FITC Conjugate (Cat. No. 8E9-FITC) was diluted in blocking buffer (1:20) and incubated over-night at 4°C. DNA was stained with Hoechst in blue. MOI = multiplicity of infection; w/o = without